





PUSH TO TAKE

contact lens case design





Contact



Chao-Mu Lin
Department of Creative Product Design
+886-975-712-049
kubomshu@hotmail.com



Jennifer Yu
Creative and Innovation Center
Tel. +886-4-2332454 ext. 6433
Fax. +886-4-23321013
jenny@asia.edu.tw
500, Luchang Rd., Wufeng, Taichung 41354, Taiwan
Patent No. TAIWAN 101203300

In the situation in which you don't have your tweezers with you when you want to remove your contact lens, you probably would use your fingers instead. Your fingers, however, might get your lens scratched.

This design allows the extraction of the lenses from the case easier by hold on the device on the sides. Using this device to lift the contact lenses from the water solution is easier. After the case is capped, the silicon rubber gels from the top and bottom are pushed together so that the case is securely closed.



Pushing Down